

# PT. KRISAN KALINDO UTAMA

Distributor Alat Kesehatan & Laboratorium Terpercaya

Email: info@krisankalindoutama.com | Telp: +62-21-826-58-722

## WEIGHING TABLE 90(L)x60(W)x80(H)



**SKU : MT-01**

**Kategori : Uncategorized**

### DESKRIPSI PRODUK

The Sartorius-style laboratory weighing table measuring 90 cm (l) × 60 cm (w) × 80 cm (h) is a specialized anti-vibration workstation. It is designed to stabilize micro, semi-micro, and analytical balances by neutralizing disturbances from air currents and structural vibrations.



Scan QR Code atau kunjungi link berikut untuk melihat detail di website:

<https://www.krisankalindoutama.com/products/weighing-table-90lx60wx80h>